

# Transcription Factor Binding To IGHM Enhancer 3 (TFE3), ELISA Kit

Cat PEK02646

---

## Species

Plant

## Full Product Name

Plant Transcription Factor Binding To IGHM Enhancer 3 (TFE3) ELISA Kit

## Product Gene Name

TFE3 elisa kit

## Molecular Weight

11,363 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Transcription factor

## NCBI Accession #

AAH26027.1

## NCBI GI #

19684176

## NCBI GeneID

7030

## NCBI Official Full Name

Transcription factor binding to IGHM enhancer 3

## NCBI Official Symbol

TFE3

## NCBI Official Synonym Full Names

transcription factor binding to IGHM enhancer 3

## NCBI Official Synonym Symbols

TFEA; RCCP2; RCCX1; bHLHe33

## NCBI Protein Information

transcription factor E3

## UniProt Gene Name

TFE3

## UniProt Synonym Gene Names

BHLHE33; bHLHe33

## UniProt Protein Name

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

---

# Transcription Factor Binding To IGHM Enhancer 3 (TFE3), ELISA Kit

Cat PEK02646

---

Transcription factor E3

## UniProt Synonym Protein Names

Class E basic helix-loop-helix protein 33; bHLHe33

## UniProt Entry Name

TFE3\_HUMAN

## UniProt Secondary Accession #

Q5JU74; Q92757; Q92758; Q99964; A8MZL6

## UniProt Related Accession #

P19532

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY