

Sterile Alpha and TIR Motif-Containing Protein 1 (SARM1), ELISA Kit

Cat PEK02649

Species

Plant

Full Product Name

Plant Sterile Alpha and TIR Motif-Containing Protein 1 (SARM1) ELISA Kit

Product Gene Name

SARM1 elisa kit

Molecular Weight

75,337 Da

Storage

Store all reagents at 2-8°C

Protein Family

Sterile alpha and TIR motif-containing protein

NCBI Accession

NP_055892.2

NCBI GI

154090976

NCBI GeneID

23098

NCBI Official Full Name

sterile alpha and TIR motif-containing protein 1

NCBI Official Symbol

SARM1

NCBI Official Synonym Full Names

sterile alpha and TIR motif containing 1

NCBI Official Synonym Symbols

SARM; SAMD2; MyD88-5

NCBI Protein Information

sterile alpha and TIR motif-containing protein 1

UniProt Gene Name

SARM1

UniProt Synonym Gene Names

KIAA0524; SAMD2; SARM; MyD88-5; SAM domain-containing protein 2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Sterile Alpha and TIR Motif-Containing Protein 1 (SARM1), ELISA Kit

Cat PEK02649

UniProt Protein Name

Sterile alpha and TIR motif-containing protein 1

UniProt Synonym Protein Names

Sterile alpha and Armadillo repeat protein; Sterile alpha motif domain-containing protein 2; MyD88-5; SAM domain-containing protein 2; Tir-1 homolog

UniProt Entry Name

SARM1_HUMAN

UniProt Secondary Accession

O60277; Q7LGG3; Q9NXY5

UniProt Related Accession

Q6SZW1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY