

Sphingomyelin Phosphodiesterase 3 (SMPD3), ELISA Kit

Cat PEK02650

Species

Plant

Full Product Name

Plant Sphingomyelin Phosphodiesterase 3 (SMPD3) ELISA Kit

Product Gene Name

SMPD3 elisa kit

Molecular Weight

70,170 Da

Storage

Store all reagents at 2-8°C

Protein Family

Sphingomyelin phosphodiesterase

NCBI Accession

NP_061137.1

NCBI GI

8923946

NCBI GeneID

55512

NCBI Official Full Name

sphingomyelin phosphodiesterase 3

NCBI Official Symbol

SMPD3

NCBI Official Synonym Full Names

sphingomyelin phosphodiesterase 3

NCBI Official Synonym Symbols

NSMASE2

NCBI Protein Information

sphingomyelin phosphodiesterase 3

UniProt Gene Name

SMPD3

UniProt Synonym Gene Names

nSMase-2; nSMase2

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Sphingomyelin Phosphodiesterase 3 (SMPD3), ELISA Kit

Cat PEK02650

Sphingomyelin phosphodiesterase 3

UniProt Synonym Protein Names

Neutral sphingomyelinase 2; nSMase-2; nSMase2; Neutral sphingomyelinase II

UniProt Entry Name

NSMA2_HUMAN

UniProt Secondary Accession

Q2M1S8; B7ZL82

UniProt Related Accession

Q9NY59

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY