

Prokineticin Receptor 2 (PKR2), ELISA Kit

Cat PEK02657

Species

Plant

Full Product Name

Plant Prokineticin Receptor 2 (PKR2) ELISA Kit

Product Gene Name

PKR2 elisa kit

Molecular Weight

43,996 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

NP_658986.1

NCBI GI

21426829

NCBI GenelD

128674

NCBI Official Full Name

prokineticin receptor 2

NCBI Official Symbol

PROKR2

NCBI Official Synonym Full Names

prokineticin receptor 2

NCBI Official Synonym Symbols

HH3; KAL3; PKR2; GPRg2; GPR73b; GPR73L1; dJ680N4.3

NCBI Protein Information

prokineticin receptor 2

UniProt Gene Name

PROKR2

UniProt Synonym Gene Names

GPR73L1; PKR2; PK-R2

UniProt Protein Name

Prokineticin receptor 2

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Prokineticin Receptor 2 (PKR2), ELISA Kit

Cat PEK02657

G-protein coupled receptor 73-like 1; G-protein coupled receptor I5E; GPR73b; GPRg2

UniProt Entry Name

PKR2_HUMAN

UniProt Secondary Accession

Q2M3C0; Q5TDY1; Q9NTT0; A5JUU1

UniProt Related Accession

Q8NFJ6

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
