

Calcium Channel, Voltage Dependent, L-Type, Alpha 1C Subunit (CACNA1C), ELISA Kit

Cat PEK02740

Species

Plant

Full Product Name

Plant Calcium Channel, Voltage Dependent, L-Type, Alpha 1C Subunit (CACNA1C) ELISA Kit

Product Gene Name

CACNA1C elisa kit

Storage

Store all reagents at 2-8°C

Protein Family

Voltage-dependent L-type calcium channel

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY