

S100 Calcium Binding Protein A5 (S100A5), ELISA Kit

Cat PEK02758

Species

Plant

Full Product Name

Plant S100 Calcium Binding Protein A5 (S100A5) ELISA Kit

Product Gene Name

S100A5 elisa kit

Molecular Weight

12,813 Da

Storage

Store all reagents at 2-8°C

Protein Family

Protein

NCBI Accession

AAH93957.1

NCBI GI

62739837

NCBI GeneID

6276

NCBI Official Full Name

S100 calcium binding protein A5

NCBI Official Symbol

S100A5

NCBI Official Synonym Full Names

S100 calcium binding protein A5

NCBI Official Synonym Symbols

S100D

NCBI Protein Information

protein S100-A5

UniProt Gene Name

S100A5

UniProt Synonym Gene Names

S100D

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

S100 Calcium Binding Protein A5 (S100A5), ELISA Kit

Cat PEK02758

Protein S100-A5

UniProt Synonym Protein Names

Protein S-100D; S100 calcium-binding protein A5

UniProt Entry Name

S10A5_HUMAN

UniProt Secondary Accession

Q52LE7; Q5RHS3

UniProt Related Accession

P33763

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY