

Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM), ELISA Kit

Cat PEK02763

Species

Plant

Full Product Name

Plant Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM) ELISA Kit

Product Gene Name

PICALM elisa kit

Molecular Weight

70,037 Da

Storage

Store all reagents at 2-8°C

Protein Family

Phosphatidylinositol-binding clathrin assembly protein

NCBI Accession

AAH48259.2

NCBI GI

57169203

NCBI GeneID

8301

NCBI Official Full Name

Phosphatidylinositol binding clathrin assembly protein

NCBI Official Symbol

PICALM

NCBI Official Synonym Full Names

phosphatidylinositol binding clathrin assembly protein

NCBI Official Synonym Symbols

LAP; CALM; CLTH

NCBI Protein Information

phosphatidylinositol-binding clathrin assembly protein

UniProt Gene Name

PICALM

UniProt Synonym Gene Names

CALM

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM), ELISA Kit

Cat PEK02763

UniProt Protein Name

Phosphatidylinositol-binding clathrin assembly protein

UniProt Synonym Protein Names

Clathrin assembly lymphoid myeloid leukemia protein

UniProt Entry Name

PICAL_HUMAN

UniProt Secondary Accession

O60700; Q4LE54; Q6GMQ6; Q86XZ9; B4DTM3; E9PN05; F8VPG7

UniProt Related Accession

Q13492

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY