

Papilin, Proteoglycan-Like Sulfated Glycoprotein (PAPLN), ELISA Kit

Cat PEK02765

Species

Plant

Full Product Name

Plant Papilin, Proteoglycan-Like Sulfated Glycoprotein (PAPLN) ELISA Kit

Product Gene Name

PAPLN elisa kit

Molecular Weight

119,209 Da

Storage

Store all reagents at 2-8°C

Protein Family

Papilin

NCBI Accession

NP_775733.3

NCBI GI

145309328

NCBI GenelD

89932

NCBI Official Full Name

papilin

NCBI Official Symbol

PAPLN

NCBI Official Synonym Full Names

papilin, proteoglycan like sulfated glycoprotein

NCBI Protein Information

papilin

UniProt Gene Name

PAPLN

UniProt Protein Name

Papilin

UniProt Entry Name

PPN_HUMAN

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Papilin, Proteoglycan-Like Sulfated Glycoprotein (PAPLN), ELISA Kit

Cat PEK02765



UniProt Secondary Accession

Q659F2; Q6UXJ4; Q6ZNM1; Q6ZUJ0; Q7Z681; Q8IVU0; B4DES8; B4DGE6

UniProt Related Accession

O95428

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY