

ISL LIM Homeobox Protein 1 (ISL1), ELISA Kit

Cat PEK02773

Species

Plant

Full Product Name

Plant ISL LIM Homeobox Protein 1 (ISL1) ELISA Kit

Product Gene Name

ISL1 elisa kit

Molecular Weight

39,036 Da

Storage

Store all reagents at 2-8°C

Protein Family

Insulin gene enhancer protein

NCBI Accession

AAH31213.1

NCBI GI

22658390

NCBI GeneID

3670

NCBI Official Full Name

ISL LIM homeobox 1

NCBI Official Symbol

ISL1

NCBI Official Synonym Full Names

ISL LIM homeobox 1

NCBI Official Synonym Symbols

Isl-1; ISLET1

NCBI Protein Information

insulin gene enhancer protein ISL-1

UniProt Gene Name

ISL1

UniProt Synonym Gene Names

Islet-1

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

ISL LIM Homeobox Protein 1 (ISL1), ELISA Kit

Cat PEK02773

Insulin gene enhancer protein ISL-1

UniProt Entry Name

ISL1_HUMAN

UniProt Secondary Accession

P20663; P47894

UniProt Related Accession

P61371

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY