

Cartilage Glycoprotein 39/Ykl40 (YKL40), ELISA Kit

Cat PEK02787

Species

Plant

Full Product Name

Plant Cartilage Glycoprotein 39/Ykl40 (YKL40) ELISA Kit

Product Gene Name

YKL40 elisa kit

Assay Type

Sandwich (Quantitative)

Detection Range

6.25ng/ml-200ng/ml.

Sensitivity

1.0ng/ml.

Intra-assay Precision

Intra-assay CV (%) is less than 15%.

Inter-assay Precision

Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean × 100].

Specificity

No significant cross-reactivity or interference between this analyte and analogues is observed.

Samples

Undiluted original Plant Tissue Homogenate and Biological Fluids samples.

Storage

Store all reagents at 2-8°C

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY