

Transcription Elongation Regulator 1-Like Protein (TCERG1L), ELISA Kit

Cat PEK02800

Species

Plant

Full Product Name

Plant Transcription Elongation Regulator 1-Like Protein (TCERG1L) ELISA Kit

Product Gene Name

TCERG1L elisa kit

Molecular Weight

65,660 Da

Storage

Store all reagents at 2-8°C

Protein Family

Transcription elongation regulator 1-like protein

NCBI Accession

NP_777597.2

NCBI GI

148727264

NCBI GenelD

256536

NCBI Official Full Name

transcription elongation regulator 1-like protein

NCBI Official Symbol

TCERG1L

NCBI Official Synonym Full Names

transcription elongation regulator 1 like

NCBI Protein Information

transcription elongation regulator 1-like protein

UniProt Gene Name

TCERG1L

UniProt Protein Name

Transcription elongation regulator 1-like protein

UniProt Entry Name

TCRGL_HUMAN

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Transcription Elongation Regulator 1-Like Protein (TCERG1L), ELISA Kit

Cat *PEK02800*

UniProt Secondary Accession

Q5VWI2; Q86XM8

UniProt Related Accession

Q5VWI1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY