

Protease Activated Receptor 2 (PAR2), ELISA Kit

Cat PEK03255

Species

Plant

Full Product Name

Plant Protease Activated Receptor 2 (PAR2) ELISA Kit

Product Gene Name

PAR2 elisa kit

Molecular Weight

44,126 Da

Storage

Store all reagents at 2-8°C

Protein Family

Transcription factor

NCBI Accession

NP_005233.3

NCBI GI

34577052

NCBI GeneID

2150

NCBI Official Full Name

proteinase-activated receptor 2

NCBI Official Symbol

F2RL1

NCBI Official Synonym Full Names

F2R like trypsin receptor 1

NCBI Official Synonym Symbols

PAR2; GPR11

NCBI Protein Information

proteinase-activated receptor 2

UniProt Gene Name

F2RL1

UniProt Synonym Gene Names

GPR11; PAR2; PAR-2

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Protease Activated Receptor 2 (PAR2), ELISA Kit

Cat PEK03255

Proteinase-activated receptor 2

UniProt Synonym Protein Names

Coagulation factor II receptor-like 1; G-protein coupled receptor 11; Thrombin receptor-like 1

UniProt Entry Name

PAR2_HUMAN

UniProt Secondary Accession

Q13317; Q13346; Q53XJ8

UniProt Related Accession

P55085

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY