

PH Domain Leucine-Rich Repeat-Containing Protein Phosphatase 1 (PHLPP), ELISA Kit

Cat PEK03937

Species

Plant

Full Product Name

Plant PH Domain Leucine-Rich Repeat-Containing Protein Phosphatase 1 (PHLPP) ELISA Kit

Product Gene Name

PHLPP elisa kit

Molecular Weight

133,565 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

NP_919431.2

NCBI GI

291219891

NCBI GeneID

23239

NCBI Official Full Name

PH domain leucine-rich repeat-containing protein phosphatase 1

NCBI Official Symbol

PHLPP1

NCBI Official Synonym Full Names

PH domain and leucine rich repeat protein phosphatase 1

NCBI Official Synonym Symbols

SCOP; PHLPP; PPM3A; PLEKHE1

NCBI Protein Information

PH domain leucine-rich repeat-containing protein phosphatase 1

UniProt Gene Name

PHLPP1

UniProt Synonym Gene Names

KIAA0606; PHLPP; PLEKHE1; SCOP; PH domain-containing family E member 1; hSCOP

UniProt Protein Name

PH domain leucine-rich repeat-containing protein phosphatase 1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

PH Domain Leucine-Rich Repeat-Containing Protein Phosphatase 1 (PHLPP), ELISA Kit

Cat PEK03937

UniProt Synonym Protein Names

Pleckstrin homology domain-containing family E member 1; PH domain-containing family E member 1; Suprachiasmatic nucleus circadian oscillatory protein; hSCOP

UniProt Entry Name

PHLP1_HUMAN

UniProt Secondary Accession

Q641Q7; Q6P4C4; Q6PJI6; Q86TN6; Q96FK2; Q9NUY1; A1A4F5

UniProt Related Accession

O60346

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY