

Lumican (LUM), ELISA Kit

Cat PEK03941

Species

Plant

Full Product Name

Plant Lumican (LUM) ELISA Kit

Product Gene Name

LUM elisa kit

Molecular Weight

38,429 Da

Storage

Store all reagents at 2-8°C

Protein Family

Lumican

NCBI Accession

NP_002336.1

NCBI GI

4505047

NCBI GeneID

4060

NCBI Official Full Name

lumican

NCBI Official Symbol

LUM

NCBI Official Synonym Full Names

lumican

NCBI Official Synonym Symbols

LDC; SLRR2D

NCBI Protein Information

lumican

UniProt Gene Name

LUM

UniProt Synonym Gene Names

LDC; SLRR2D; KSPG lumican

UniProt Protein Name

Lumican

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Lumican (LUM), ELISA Kit

Cat PEK03941

UniProt Synonym Protein Names

Keratan sulfate proteoglycan lumican; KSPG lumican

UniProt Entry Name

LUM_HUMAN

UniProt Secondary Accession

Q96QM7; B2R6R5

UniProt Related Accession

P51884

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY