

# Catechol (CT), ELISA Kit

Cat PEK03969

---

## Species

Plant

## Full Product Name

Plant Catechol (CT) ELISA Kit

## Product Gene Name

CT elisa kit

## Assay Type

Sandwich

## Detection Range

31.2ng/ml-1000ng/ml.

## Sensitivity

5.0ng/ml.

## Intra-assay Precision

Intra-assay CV (%) is less than 15%.

## Inter-assay Precision

Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean × 100]

## Storage

Store all reagents at 2-8°C

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**