

C-C Chemokine Receptor Type 8 (CCR8), ELISA Kit

Cat PEK03971

Species

Plant

Full Product Name

Plant C-C Chemokine Receptor Type 8 (CCR8) ELISA Kit

Product Gene Name

CCR8 elisa kit

Molecular Weight

31,622 Da

Storage

Store all reagents at 2-8°C

Protein Family

C-C chemokine receptor

NCBI Accession

NP_005192.1

NCBI GI

4885121

NCBI GeneID

1237

NCBI Official Full Name

C-C chemokine receptor type 8

NCBI Official Symbol

CCR8

NCBI Official Synonym Full Names

C-C motif chemokine receptor 8

NCBI Official Synonym Symbols

CY6; TER1; CCR-8; CKRL1; CDw198; CMKBR8; GPRCY6; CMKBRL2; CC-CKR-8

NCBI Protein Information

C-C chemokine receptor type 8

UniProt Gene Name

CCR8

UniProt Synonym Gene Names

CKRL1; CMKBR8; CMKBRL2; C-C CKR-8; CC-CKR-8; CCR-8; CKR-L1; GPRCY6

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

C-C Chemokine Receptor Type 8 (CCR8), ELISA Kit

Cat PEK03971

C-C chemokine receptor type 8

UniProt Synonym Protein Names

CC chemokine receptor CHEMR1; CMKBRL2; Chemokine receptor-like 1; CKR-L1; GPR-CY6; GPRCY6; TER1;
CD_antigen: CDw198

UniProt Entry Name

CCR8_HUMAN

UniProt Secondary Accession

Q3KNQ8; Q3KNR3; Q9BYX5; B2RC64

UniProt Related Accession

P51685

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY