

# Teashirt Homolog 3 (TSHZ3), ELISA Kit

Cat PEK03985

---

## Species

Plant

## Full Product Name

Plant Teashirt Homolog 3 (TSHZ3) ELISA Kit

## Product Gene Name

TSHZ3 elisa kit

## Molecular Weight

118,566 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Teashirt

## NCBI Accession #

NP\_065907.2

## NCBI GI #

127138957

## NCBI GeneID

57616

## NCBI Official Full Name

teashirt homolog 3

## NCBI Official Symbol

TSHZ3

## NCBI Official Synonym Full Names

teashirt zinc finger homeobox 3

## NCBI Official Synonym Symbols

TSH3; ZNF537

## NCBI Protein Information

teashirt homolog 3

## UniProt Gene Name

TSHZ3

## UniProt Synonym Gene Names

KIAA1474; TSH3; ZNF537

## UniProt Protein Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Teashirt Homolog 3 (TSHZ3), ELISA Kit

Cat PEK03985

---

Teashirt homolog 3

## UniProt Synonym Protein Names

Zinc finger protein 537

## UniProt Entry Name

TSH3\_HUMAN

## UniProt Secondary Accession #

Q9H0G6; Q9P254

## UniProt Related Accession #

Q63HK5

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**