

Protein Tyrosine Phosphatase Receptor Type H (PTPRH), ELISA Kit

Cat PEK03999

Species

Plant

Full Product Name

Plant Protein Tyrosine Phosphatase Receptor Type H (PTPRH) ELISA Kit

Product Gene Name

PTPRH elisa kit

Molecular Weight

103,412 Da

Storage

Store all reagents at 2-8°C

Protein Family

Receptor-type tyrosine-protein phosphatase

NCBI Accession

AAI11717.1

NCBI GI

84627436

NCBI GeneID

5794

NCBI Official Full Name

Protein tyrosine phosphatase, receptor type, H

NCBI Official Symbol

PTPRH

NCBI Official Synonym Full Names

protein tyrosine phosphatase, receptor type H

NCBI Official Synonym Symbols

SAP1; R-PTP-H

NCBI Protein Information

receptor-type tyrosine-protein phosphatase H

UniProt Gene Name

PTPRH

UniProt Synonym Gene Names

SAP1; R-PTP-H; SAP-1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Protein Tyrosine Phosphatase Receptor Type H (PTPRH), ELISA Kit

Cat PEK03999

UniProt Protein Name

Receptor-type tyrosine-protein phosphatase H

UniProt Synonym Protein Names

Stomach cancer-associated protein tyrosine phosphatase 1; SAP-1; Transmembrane-type protein-tyrosine phosphatase type H

UniProt Entry Name

PTPRH_HUMAN

UniProt Secondary Accession

Q15426; Q2NKN9; Q2NKP0; C9JCH2

UniProt Related Accession

Q9HD43

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY