

Leukotriene A4 Hydrolase (LTA4H), ELISA Kit

Cat PEK04011

Species

Plant

Full Product Name

Plant Leukotriene A4 Hydrolase (LTA4H) ELISA Kit

Product Gene Name

LTA4H elisa kit

Molecular Weight

66,852 Da

Storage

Store all reagents at 2-8°C

Protein Family

Leukotriene A-4 hydrolase

NCBI Accession

AAA89077.1

NCBI GI

976396

NCBI Official Full Name

leukotriene A4 hydrolase

UniProt Gene Name

LTA4H

UniProt Synonym Gene Names

LTA4; LTA-4 hydrolase

UniProt Protein Name

Leukotriene A-4 hydrolase

UniProt Synonym Protein Names

Leukotriene A(4) hydrolase

UniProt Entry Name

LKHA4_HUMAN

UniProt Secondary Accession

Q6IAT6; Q9UCT7; B4DNQ9; F8VV40

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY