

Interleukin 4 Induced Protein 1 (IL4I1), ELISA Kit

Cat PEK04013

Species

Plant

Full Product Name

Plant Interleukin 4 Induced Protein 1 (IL4I1) ELISA Kit

Product Gene Name

IL4I1 elisa kit

Molecular Weight

65,328 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

CAI54291.1

NCBI GI

71060509

NCBI GenelD

259307

NCBI Official Full Name

interleukin 4 induced protein 1

NCBI Official Symbol

IL4I1

NCBI Official Synonym Full Names

interleukin 4 induced 1

NCBI Official Synonym Symbols

LAO; FIG1; LAAO

NCBI Protein Information

L-amino-acid oxidase

UniProt Gene Name

IL4I1

UniProt Synonym Gene Names

FIG1; LAAO; LAO; IL4-induced protein 1; hFIG1

UniProt Protein Name

L-amino-acid oxidase

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Interleukin 4 Induced Protein 1 (IL4I1), ELISA Kit

Cat PEK04013

Interleukin-4-induced protein 1; IL4-induced protein 1; Protein Fig-1; hFIG1

UniProt Entry Name

OXLA_HUMAN

UniProt Primary Accession

Q96RQ9

UniProt Secondary Accession

Q1WMJ3; Q4GZN1; Q4GZN2; Q6P2Q3; Q8TEM5; Q96RQ8

UniProt Related Accession

Q96RQ9

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY