

BAG Family Molecular Chaperone Regulator 2 (BAG2), ELISA Kit

Cat PEK04037

Species

Plant

Full Product Name

Plant BAG Family Molecular Chaperone Regulator 2 (BAG2) ELISA Kit

Product Gene Name

BAG2 elisa kit

Molecular Weight

20,189 Da

Storage

Store all reagents at 2-8°C

Protein Family

BAG family molecular chaperone regulator

NCBI Accession

NP_004273.1

NCBI GI

4757834

NCBI GeneID

9532

NCBI Official Full Name

BAG family molecular chaperone regulator 2

NCBI Official Symbol

BAG2

NCBI Official Synonym Full Names

BCL2 associated athanogene 2

NCBI Official Synonym Symbols

BAG-2; dJ417I1.2

NCBI Protein Information

BAG family molecular chaperone regulator 2

UniProt Gene Name

BAG2

UniProt Synonym Gene Names

BAG-2

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

BAG Family Molecular Chaperone Regulator 2 (BAG2), ELISA Kit

Cat PEK04037

BAG family molecular chaperone regulator 2

UniProt Synonym Protein Names

Bcl-2-associated athanogene 2

UniProt Entry Name

BAG2_HUMAN

UniProt Secondary Accession

Q08AS9; Q6FID0; B4DXE2

UniProt Related Accession

O95816

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
