

Uncoupling Protein 3 (UCP3), ELISA Kit

Cat PEK04185

Species

Plant

Full Product Name

Plant Uncoupling Protein 3 (UCP3) ELISA Kit

Product Gene Name

UCP3 elisa kit

Molecular Weight

22,874 Da

Storage

Store all reagents at 2-8°C

Protein Family

UBA domain-containing protein

NCBI Accession

AAC18822.1

NCBI GI

3176760

NCBI Official Full Name

uncoupling protein 3

UniProt Gene Name

UCP3

UniProt Synonym Gene Names

SLC25A9; UCP 3

UniProt Protein Name

Mitochondrial uncoupling protein 3

UniProt Synonym Protein Names

Solute carrier family 25 member 9

UniProt Entry Name

UCP3_HUMAN

UniProt Secondary Accession

O60475; Q96HL3

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY