

Plasminogen Activator Inhibitor 2 (PAI2), ELISA Kit

Cat PEK04280

Species

Plant

Full Product Name

Plant Plasminogen Activator Inhibitor 2 (PAI2) ELISA Kit

Product Gene Name

PAI2 elisa kit

Molecular Weight

46,596 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

NP_001137290.1

NCBI GI

219689110

NCBI GeneID

5055

NCBI Official Full Name

plasminogen activator inhibitor 2

NCBI Official Symbol

SERPINB2

NCBI Official Synonym Full Names

serpin family B member 2

NCBI Official Synonym Symbols

PAI; PAI2; PAI-2; PLANH2; HsT1201

NCBI Protein Information

plasminogen activator inhibitor 2

UniProt Gene Name

SERPINB2

UniProt Synonym Gene Names

PAI2; PLANH2; PAI-2

UniProt Protein Name

Plasminogen activator inhibitor 2

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Plasminogen Activator Inhibitor 2 (PAI2), ELISA Kit

Cat PEK04280

Monocyte Arg-serpin; Placental plasminogen activator inhibitor; Serpin B2; Urokinase inhibitor

UniProt Entry Name

PAI2_HUMAN

UniProt Secondary Accession

Q96E96

UniProt Related Accession

P05120

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY