

Low Density Lipoprotein (LDL), ELISA Kit

Cat PEK04297

Species

Plant

Full Product Name

Plant Low Density Lipoprotein (LDL) ELISA Kit

Product Gene Name

LDL elisa kit

Molecular Weight

82,255 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

ADM85894.1

NCBI GI

306411097

NCBI GenID

3949

NCBI Official Full Name

low density lipoprotein, partial

NCBI Official Symbol

LDLR

NCBI Official Synonym Full Names

low density lipoprotein receptor

NCBI Official Synonym Symbols

FH; FHC; LDLCQ2

NCBI Protein Information

low-density lipoprotein receptor

UniProt Gene Name

LDLR

UniProt Synonym Gene Names

LDL receptor

UniProt Protein Name

Low-density lipoprotein receptor

UniProt Entry Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Low Density Lipoprotein (LDL), ELISA Kit

Cat PEK04297

LDLR_HUMAN

UniProt Secondary Accession

Q53ZD9; Q59FQ1; Q9UDH7; B4DII3; B4DJZ8; B4DR00; B4DTQ3; C0JYY8; H0YLU8; H0YNT7

UniProt Related Accession

P01130

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY