

Legumain (LGMN), ELISA Kit

Cat PEK04300

Species

Plant

Full Product Name

Plant Legumain (LGMN) ELISA Kit

Product Gene Name

LGMN elisa kit

Molecular Weight

42,008 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

NP_005597.3

NCBI GI

56682962

NCBI GeneID

5641

NCBI Official Full Name

legumain preproprotein

NCBI Official Symbol

LGMN

NCBI Official Synonym Full Names

legumain

NCBI Official Synonym Symbols

AEP; LGMN1; PRSC1

NCBI Protein Information

legumain

UniProt Gene Name

LGMN

UniProt Synonym Gene Names

PRSC1

UniProt Protein Name

Legumain

UniProt Synonym Protein Names

Asparaginyl endopeptidase; Protease, cysteine 1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Legumain (LGMN), ELISA Kit

Cat PEK04300

UniProt Entry Name

LGMN_HUMAN

UniProt Secondary Accession

O00123; Q86TV2; Q86TV3; Q9BTY1

UniProt Related Accession

Q99538

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY