

Guanine Nucleotide-Binding Protein Subunit Beta-4 (GNB4), ELISA Kit

Cat PEK04312

Species

Plant

Full Product Name

Plant Guanine Nucleotide-Binding Protein Subunit Beta-4 (GNB4) ELISA Kit

Product Gene Name

GNB4 elisa kit

Molecular Weight

37,567 Da

Storage

Store all reagents at 2-8°C

Protein Family

Guanine nucleotide-binding protein

NCBI Accession

NP_067642.1

NCBI GI

11055998

NCBI GenelD

59345

NCBI Official Full Name

guanine nucleotide-binding protein subunit beta-4

NCBI Official Symbol

GNB4

NCBI Official Synonym Full Names

G protein subunit beta 4

NCBI Official Synonym Symbols

CMTD1F

NCBI Protein Information

guanine nucleotide-binding protein subunit beta-4

UniProt Gene Name

GNB4

UniProt Protein Name

Guanine nucleotide-binding protein subunit beta-4

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Guanine Nucleotide-Binding Protein Subunit Beta-4 (GNB4), ELISA Kit

Cat PEK04312

UniProt Synonym Protein Names

Transducin beta chain 4

UniProt Entry Name

GBB4_HUMAN

UniProt Secondary Accession

B3KMH5; D3DNR8

UniProt Related Accession

Q9HAV0

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
