

# Endothelin 1 (EDN1), ELISA Kit

Cat PEK04320

## Species

Plant

## Full Product Name

Plant Endothelin 1 (EDN1) ELISA Kit

## Product Gene Name

EDN1 elisa kit

## Assay Type

Quantitative Sandwich

## Detection Range

6.25pg/ml-200pg/mL

## Sensitivity

1.0pg/mL

## Intra-assay Precision

Intra-assay CV (%) is less than 15%.

## Inter-assay Precision

Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean × 100].

## Storage

Store all reagents at 2-8°C

## Protein Family

Endothelin

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**