

# Sphingosine 1 Phosphate Receptor 5 (S1PR5), ELISA Kit

Cat PEK04408

---

## Species

Plant

## Full Product Name

Plant Sphingosine 1 Phosphate Receptor 5 (S1PR5) ELISA Kit

## Product Gene Name

S1PR5 elisa kit

## Molecular Weight

32,279 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Sphingosine 1-phosphate receptor

## NCBI Accession #

NP\_001159687.1

## NCBI GI #

261278306

## NCBI GeneID

53637

## NCBI Official Full Name

sphingosine 1-phosphate receptor 5

## NCBI Official Symbol

S1PR5

## NCBI Official Synonym Full Names

sphingosine-1-phosphate receptor 5

## NCBI Official Synonym Symbols

EDG8; S1P5; Edg-8; SPPR-1; SPPR-2

## NCBI Protein Information

sphingosine 1-phosphate receptor 5

## UniProt Gene Name

S1PR5

## UniProt Synonym Gene Names

EDG8; S1P receptor 5; S1P5; S1P receptor Edg-8

## UniProt Protein Name

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

---

# Sphingosine 1 Phosphate Receptor 5 (S1PR5), ELISA Kit

Cat PEK04408

---

Sphingosine 1-phosphate receptor 5

## UniProt Synonym Protein Names

Endothelial differentiation G-protein-coupled receptor 8; Sphingosine 1-phosphate receptor Edg-8; S1P receptor Edg-8

## UniProt Entry Name

S1PR5\_HUMAN

## UniProt Secondary Accession #

Q6NW11

## UniProt Related Accession #

Q9H228

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**