

Semaphorin 3C (SEMA3C), ELISA Kit

Cat PEK04411

Species

Plant

Full Product Name

Plant Semaphorin 3C (SEMA3C) ELISA Kit

Product Gene Name

SEMA3C elisa kit

Molecular Weight

26,255 Da

Storage

Store all reagents at 2-8°C

Protein Family

Semaphorin

NCBI Accession

NP_006370.1

NCBI GI

5454048

NCBI GeneID

10512

NCBI Official Full Name

semaphorin-3C

NCBI Official Symbol

SEMA3C

NCBI Official Synonym Full Names

semaphorin 3C

NCBI Official Synonym Symbols

SemE; SEMAE

NCBI Protein Information

semaphorin-3C

UniProt Gene Name

SEMA3C

UniProt Synonym Gene Names

SEMAE; Sema E

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Semaphorin 3C (SEMA3C), ELISA Kit

Cat PEK04411

Semaphorin-3C

UniProt Synonym Protein Names

Semaphorin-E; Sema E

UniProt Entry Name

SEM3C_HUMAN

UniProt Secondary Accession

B4DRL8

UniProt Related Accession

Q99985

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY