

Protein Tyrosine Phosphatase Receptor Type O (PTPRO), ELISA Kit

Cat PEK04413

Species

Plant

Full Product Name

Plant Protein Tyrosine Phosphatase Receptor Type O (PTPRO) ELISA Kit

Product Gene Name

PTPRO elisa kit

Molecular Weight

67,672 Da

Storage

Store all reagents at 2-8°C

Protein Family

Receptor-type tyrosine-protein phosphatase

NCBI Accession

AAI26204.1

NCBI GI

116496981

NCBI GeneID

5800

NCBI Official Full Name

Protein tyrosine phosphatase, receptor type, O

NCBI Official Symbol

PTPRO

NCBI Official Synonym Full Names

protein tyrosine phosphatase, receptor type O

NCBI Official Synonym Symbols

NPHS6; PTPU2; GLEPP1; PTP-OC; PTP-U2; PTPROT; R-PTP-O

NCBI Protein Information

receptor-type tyrosine-protein phosphatase O

UniProt Gene Name

PTPRO

UniProt Synonym Gene Names

GLEPP1; PTPU2; R-PTP-O; PTP-U2; PTPase U2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Protein Tyrosine Phosphatase Receptor Type O (PTPRO), ELISA Kit

Cat PEK04413

UniProt Protein Name

Receptor-type tyrosine-protein phosphatase O

UniProt Synonym Protein Names

Glomerular epithelial protein 1; Protein tyrosine phosphatase U2; PTP-U2; PTPase U2

UniProt Entry Name

PTPRO_HUMAN

UniProt Secondary Accession

Q13101; Q8IYG3; Q9UBF0; Q9UBT5; A0AV39

UniProt Related Accession

Q16827

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY