

Dolichol Kinase (DOLK), ELISA Kit

Cat PEK04449

Species

Plant

Full Product Name

Plant Dolichol Kinase (DOLK) ELISA Kit

Product Gene Name

DOLK elisa kit

Molecular Weight

59,268 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

NP_055723.1

NCBI GI

7662482

NCBI GeneID

22845

NCBI Official Full Name

dolichol kinase

NCBI Official Symbol

DOLK

NCBI Official Synonym Full Names

dolichol kinase

NCBI Official Synonym Symbols

DK; DK1; CDG1M; SEC59; TMEM15

NCBI Protein Information

dolichol kinase

UniProt Gene Name

DOLK

UniProt Synonym Gene Names

KIAA1094; TMEM15

UniProt Protein Name

Dolichol kinase

UniProt Synonym Protein Names

Transmembrane protein 15

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Dolichol Kinase (DOLK), ELISA Kit

Cat PEK04449



UniProt Entry Name

DOLK_HUMAN

UniProt Secondary Accession

Q5SRE6

UniProt Related Accession

Q9UPQ8

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY