

# Carbonic Anhydrase III, Muscle Specific (CA3), ELISA Kit

Cat PEK04619

---

## Species

Plant

## Full Product Name

Plant Carbonic Anhydrase III, Muscle Specific (CA3) ELISA Kit

## Product Gene Name

CA3 elisa kit

## Molecular Weight

29,557 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Carbonic anhydrase

## NCBI Accession #

NP\_005172.1

## NCBI GI #

4885099

## NCBI GeneID

761

## NCBI Official Full Name

carbonic anhydrase 3

## NCBI Official Symbol

CA3

## NCBI Official Synonym Full Names

carbonic anhydrase 3

## NCBI Official Synonym Symbols

Car3; CAIII

## NCBI Protein Information

carbonic anhydrase 3

## UniProt Gene Name

CA3

## UniProt Synonym Gene Names

CA-III

## UniProt Protein Name

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

# Carbonic Anhydrase III, Muscle Specific (CA3), ELISA Kit

Cat PEK04619

---

Carbonic anhydrase 3

## UniProt Synonym Protein Names

Carbonate dehydratase III; Carbonic anhydrase III; CA-III

## UniProt Entry Name

CAH3\_HUMAN

## UniProt Secondary Accession #

O60842; B2R867; B3KUC8

## UniProt Related Accession #

P07451

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**