

# Dynein Heavy Chain 5, Axonemal (DNAH5), ELISA Kit

Cat PEK04671

---

## Species

Plant

## Full Product Name

Plant Dynein Heavy Chain 5, Axonemal (DNAH5) ELISA Kit

## Product Gene Name

DNAH5 elisa kit

## Molecular Weight

529,021 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Dynein

## NCBI Accession #

NP\_001360.1

## NCBI GI #

19115954

## NCBI GeneID

1767

## NCBI Official Full Name

dynein heavy chain 5, axonemal

## NCBI Official Symbol

DNAH5

## NCBI Official Synonym Full Names

dynein axonemal heavy chain 5

## NCBI Official Synonym Symbols

HL1; PCD; CILD3; KTGNR; DNAHC5

## NCBI Protein Information

dynein heavy chain 5, axonemal

## UniProt Gene Name

DNAH5

## UniProt Synonym Gene Names

DNAHC5; HL1; KIAA1603

## UniProt Protein Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Dynein Heavy Chain 5, Axonemal (DNAH5), ELISA Kit

Cat PEK04671

---

Dynein heavy chain 5, axonemal

## UniProt Synonym Protein Names

Axonemal beta dynein heavy chain 5; Ciliary dynein heavy chain 5

## UniProt Entry Name

DYH5\_HUMAN

## UniProt Secondary Accession #

Q92860; Q96L74; Q9H5S7; Q9HCG9

## UniProt Related Accession #

Q8TE73

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**