

# Dynamin-Like 120kDa Protein, Mitochondrial (OPA1), ELISA Kit

Cat PEK04672

---

## Species

Plant

## Full Product Name

Plant Dynamin-Like 120kDa Protein, Mitochondrial (OPA1) ELISA Kit

## Product Gene Name

OPA1 elisa kit

## Storage

Store all reagents at 2-8°C

## Protein Family

Dynamin-like 120 kDa protein

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**