

Stathmin 1 (STMN1), ELISA Kit

Cat PEK04701

Species

Plant

Full Product Name

Plant Stathmin 1 (STMN1) ELISA Kit

Product Gene Name

STMN1 elisa kit

Molecular Weight

19,824 Da

Storage

Store all reagents at 2-8°C

Protein Family

Stathmin

NCBI Accession

BAG70262.1

NCBI GI

197692597

NCBI GeneID

3925

NCBI Official Full Name

stathmin 1, partial

NCBI Official Symbol

STMN1

NCBI Official Synonym Full Names

stathmin 1

NCBI Official Synonym Symbols

Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215

NCBI Protein Information

stathmin

UniProt Gene Name

STMN1

UniProt Synonym Gene Names

C1orf215; LAP18; OP18; Op18; pp19

UniProt Protein Name

Stathmin

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Stathmin 1 (STMN1), ELISA Kit

Cat PEK04701

UniProt Synonym Protein Names

Leukemia-associated phosphoprotein p18; Metablastin; Oncoprotein 18; Op18; Phosphoprotein p19; pp19; Prosolin; Protein Pr22; pp17

UniProt Entry Name

STMN1_HUMAN

UniProt Secondary Accession

A2A2D1; B2R4E7; B7Z8N4; D3DPJ5

UniProt Related Accession

P16949

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY