

Twist-Related Protein 2 (TWIST2), ELISA Kit

Cat PEK05103

Species

Plant

Full Product Name

Plant Twist-Related Protein 2 (TWIST2) ELISA Kit

Product Gene Name

TWIST2 elisa kit

Molecular Weight

18,124 Da

Storage

Store all reagents at 2-8°C

Protein Family

Twist-related protein

NCBI Accession

NP_001258822.1

NCBI GI

429325228

NCBI GeneID

117581

NCBI Official Full Name

twist-related protein 2

NCBI Official Symbol

TWIST2

NCBI Official Synonym Full Names

twist family bHLH transcription factor 2

NCBI Official Synonym Symbols

AMS; FFDD3; BBRSAY; DERMO1; SETLSS; bHLHa39

NCBI Protein Information

twist-related protein 2

UniProt Gene Name

TWIST2

UniProt Synonym Gene Names

BHLHA39; DERMO1; bHLHa39; Dermo-1

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Twist-Related Protein 2 (TWIST2), ELISA Kit

Cat PEK05103

Twist-related protein 2

UniProt Synonym Protein Names

Class A basic helix-loop-helix protein 39; bHLHa39; Dermis-expressed protein 1; Dermo-1

UniProt Entry Name

TWST2_HUMAN

UniProt Secondary Accession

Q3SYL6

UniProt Related Accession

Q8WVJ9

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY