

Transferrin (TF), ELISA Kit

Cat PEK05106

Species

Plant

Full Product Name

Plant Transferrin (TF) ELISA Kit

Product Gene Name

TF elisa kit

Molecular Weight

77,064 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

AAA61140.1

NCBI GI

339453

NCBI GeneID

7018

NCBI Official Full Name

transferrin

NCBI Official Symbol

TF

NCBI Official Synonym Full Names

transferrin

NCBI Official Synonym Symbols

TFQTL1; PRO1557; PRO2086; HEL-S-71p

NCBI Protein Information

serotransferrin

UniProt Gene Name

TF

UniProt Synonym Gene Names

Transferrin

UniProt Protein Name

Serotransferrin

UniProt Synonym Protein Names

Beta-1 metal-binding globulin; Siderophilin

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Transferrin (TF), ELISA Kit

Cat PEK05106

UniProt Entry Name

TRFE_HUMAN

UniProt Secondary Accession

O43890; Q1HBA5; Q9NQB8; Q9UHV0

UniProt Related Accession

P02787

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY