

Talin 2 (TLN2), ELISA Kit

Cat PEK05109

Species

Plant

Full Product Name

Plant Talin 2 (TLN2) ELISA Kit

Product Gene Name

TLN2 elisa kit

Molecular Weight

271,613 Da

Storage

Store all reagents at 2-8°C

Protein Family

Talin

NCBI Accession

AAM73764.1

NCBI GI

21666571

NCBI GeneID

83660

NCBI Official Full Name

talin 2

NCBI Official Symbol

TLN2

NCBI Official Synonym Full Names

talin 2

NCBI Official Synonym Symbols

ILWEQ

NCBI Protein Information

talin-2

UniProt Gene Name

TLN2

UniProt Synonym Gene Names

KIAA0320

UniProt Protein Name

Talin-2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Talin 2 (TLN2), ELISA Kit

Cat PEK05109

UniProt Entry Name

TLN2_HUMAN

UniProt Secondary Accession

A6NLB8

UniProt Related Accession

Q9Y4G6

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY