

TU3A Protein (TU3A), ELISA Kit

Cat PEK05110

Species

Plant

Full Product Name

Plant TU3A Protein (TU3A) ELISA Kit

Product Gene Name

TU3A elisa kit

Molecular Weight

20,872 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

BAA83072.1

NCBI GI

5689225

NCBI Official Full Name

TU3A

UniProt Gene Name

FAM107A

UniProt Synonym Gene Names

DRR1; TU3A

UniProt Protein Name

Protein FAM107A

UniProt Synonym Protein Names

Down-regulated in renal cell carcinoma 1; Protein TU3A

UniProt Entry Name

F107A_HUMAN

UniProt Secondary Accession

Q96NH4; B3KNQ4; B7ZAY5; J3KR61

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY