

# Ribonuclease Inhibitor (RI), ELISA Kit

Cat PEK05115

---

## Species

Plant

## Full Product Name

Plant Ribonuclease Inhibitor (RI) ELISA Kit

## Product Gene Name

RI elisa kit

## Molecular Weight

49,973 Da

## Storage

Store all reagents at 2-8°C

## NCBI Accession #

NP\_976319.1

## NCBI GI #

42822874

## NCBI GeneID

6050

## NCBI Official Full Name

ribonuclease inhibitor

## NCBI Official Symbol

RNH1

## NCBI Official Synonym Full Names

ribonuclease/angiogenin inhibitor 1

## NCBI Official Synonym Symbols

RAI; RNH

## NCBI Protein Information

ribonuclease inhibitor

## UniProt Gene Name

RNH1

## UniProt Synonym Gene Names

PRI; RNH; Placental RNase inhibitor; RAI

## UniProt Protein Name

Ribonuclease inhibitor

## UniProt Synonym Protein Names

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Ribonuclease Inhibitor (RI), ELISA Kit

Cat PEK05115

---

Placental ribonuclease inhibitor; Placental RNase inhibitor; Ribonuclease/angiogenin inhibitor 1; RAI

## UniProt Entry Name

RINI\_HUMAN

## UniProt Secondary Accession #

Q8IZK8; Q96FD7; Q9BQ80; Q9UDK6

## UniProt Related Accession #

P13489

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**