

Latent Membrane Protein-1 (LMP-1), ELISA Kit

Cat PEK05268

Species

Plant

Full Product Name

Plant Latent Membrane Protein-1 (LMP-1) ELISA Kit

Product Gene Name

LMP-1 elisa kit

Molecular Weight

5,443 Da

Storage

Store all reagents at 2-8°C

Protein Family

LAMP family protein

NCBI Accession

AAF22681.1

NCBI GI

6653380

NCBI Official Full Name

latent membrane protein-1, partial

UniProt Gene Name

LMP-1

UniProt Protein Name

Latent membrane protein-1

UniProt Synonym Protein Names

Latent membrane protein-1

UniProt Entry Name

Q9QAM5_EBVG

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY