

Hedgehog Interacting Protein (HIP), ELISA Kit

Cat PEK05273

Species

Plant

Full Product Name

Plant Hedgehog Interacting Protein (HIP) ELISA Kit

Product Gene Name

HIP elisa kit

Molecular Weight

36,373 Da

Storage

Store all reagents at 2-8°C

Protein Family

HSP-interacting protein

NCBI Accession

NP_071920.1

NCBI GI

20143973

NCBI GeneID

64399

NCBI Official Full Name

hedgehog-interacting protein

NCBI Official Symbol

HHIP

NCBI Official Synonym Full Names

hedgehog interacting protein

NCBI Official Synonym Symbols

HIP

NCBI Protein Information

hedgehog-interacting protein

UniProt Gene Name

HHIP

UniProt Synonym Gene Names

HIP; HHIP; HIP

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Hedgehog Interacting Protein (HIP), ELISA Kit

Cat PEK05273

Hedgehog-interacting protein

UniProt Entry Name

HHIP_HUMAN

UniProt Secondary Accession

Q6PK09; Q8NCI7; Q9B XK3; Q9H1J4; Q9H7E7

UniProt Related Accession

Q96QV1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY