

D-Dopachrome Tautomerase (DDT), ELISA Kit

Cat PEK05281

Species

Plant

Full Product Name

Plant D-Dopachrome Tautomerase (DDT) ELISA Kit

Product Gene Name

DDT elisa kit

Molecular Weight

11,323 Da

Storage

Store all reagents at 2-8°C

Protein Family

D-dopachrome decarboxylase

NCBI Accession

AAC77468.1

NCBI GI

2352915

NCBI Official Full Name

D-dopachrome tautomerase

UniProt Gene Name

DDT

UniProt Protein Name

D-dopachrome decarboxylase

UniProt Synonym Protein Names

D-dopachrome tautomerase; Phenylpyruvate tautomerase II

UniProt Entry Name

DOPD_HUMAN

UniProt Secondary Accession

O00774; O60787; Q13534; B7Z522

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY