

Chemokine C-C-Motif Receptor 1 (CCR1), ELISA Kit

Cat PEK05287

Species

Plant

Full Product Name

Plant Chemokine C-C-Motif Receptor 1 (CCR1) ELISA Kit

Product Gene Name

CCR1 elisa kit

Molecular Weight

41,173 Da

Storage

Store all reagents at 2-8°C

Protein Family

C-C chemokine receptor

NCBI Accession

NP_001286.1

NCBI GI

4502631

NCBI GeneID

1230

NCBI Official Full Name

C-C chemokine receptor type 1

NCBI Official Symbol

CCR1

NCBI Official Synonym Full Names

C-C motif chemokine receptor 1

NCBI Official Synonym Symbols

CKR1; CD191; CKR-1; HM145; CMKBR1; MIP1aR; SCYAR1

NCBI Protein Information

C-C chemokine receptor type 1

UniProt Gene Name

CCR1

UniProt Synonym Gene Names

CMKBR1; CMKR1; SCYAR1; C-C CKR-1; CC-CKR-1; CCR-1; CCR1; MIP-1alpha-R

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Chemokine C-C-Motif Receptor 1 (CCR1), ELISA Kit

Cat PEK05287

C-C chemokine receptor type 1

UniProt Synonym Protein Names

HM145; LD78 receptor; Macrophage inflammatory protein 1-alpha receptor; MIP-1alpha-R; RANTES-R;
CD_antigen: CD191

UniProt Entry Name

CCR1_HUMAN

UniProt Secondary Accession

Q86VA9

UniProt Related Accession

P32246

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY