

Chemokine C-C-Motif Ligand 1 (CCL1), ELISA Kit

Cat PEK05288

Species

Plant

Full Product Name

Plant Chemokine C-C-Motif Ligand 1 (CCL1) ELISA Kit

Product Gene Name

CCL1 elisa kit

Molecular Weight

10,992 Da

Storage

Store all reagents at 2-8°C

Protein Family

Cyclin

NCBI Accession

NP_002972.1

NCBI GI

4506833

NCBI GeneID

6346

NCBI Official Full Name

C-C motif chemokine 1

NCBI Official Symbol

CCL1

NCBI Official Synonym Full Names

C-C motif chemokine ligand 1

NCBI Official Synonym Symbols

P500; SISE; TCA3; I-309; SCYA1

NCBI Protein Information

C-C motif chemokine 1

UniProt Gene Name

CCL1

UniProt Synonym Gene Names

SCYA1

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Chemokine C-C-Motif Ligand 1 (CCL1), ELISA Kit

Cat PEK05288

C-C motif chemokine 1

UniProt Synonym Protein Names

Small-inducible cytokine A1; T lymphocyte-secreted protein I-309

UniProt Entry Name

CCL1_HUMAN

UniProt Secondary Accession

Q2M309; B2R5G9

UniProt Related Accession

P22362

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY