

Autoimmune Regulator (AIRE), ELISA Kit

Cat PEK05290

Species

Plant

Full Product Name

Plant Autoimmune Regulator (AIRE) ELISA Kit

Product Gene Name

AIRE elisa kit

Molecular Weight

36,501 Da

Storage

Store all reagents at 2-8°C

Protein Family

Autoimmune regulator

NCBI Accession

NP_000374.1

NCBI GI

4557291

NCBI GeneID

326

NCBI Official Full Name

autoimmune regulator

NCBI Official Symbol

AIRE

NCBI Official Synonym Full Names

autoimmune regulator

NCBI Official Synonym Symbols

APS1; APSI; PGA1; AIRE1; APECED

NCBI Protein Information

autoimmune regulator

UniProt Gene Name

AIRE

UniProt Synonym Gene Names

APECED; APECED protein

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Autoimmune Regulator (AIRE), ELISA Kit

Cat PEK05290

Autoimmune regulator

UniProt Synonym Protein Names

Autoimmune polyendocrinopathy candidiasis ectodermal dystrophy protein; APECED protein

UniProt Entry Name

AIRE_HUMAN

UniProt Secondary Accession

O43922; O43932; O75745; B2RP50

UniProt Related Accession

O43918

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY