

Sialidase-3 (NEU3), ELISA Kit

Cat PEK05299

Species

Plant

Full Product Name

Plant Sialidase-3 (NEU3) ELISA Kit

Product Gene Name

NEU3 elisa kit

Molecular Weight

51,675 Da

Storage

Store all reagents at 2-8°C

Protein Family

Sialidase

NCBI Accession

NP_006647.3

NCBI GI

117190519

NCBI GeneID

10825

NCBI Official Full Name

sialidase-3

NCBI Official Symbol

NEU3

NCBI Official Synonym Full Names

neuraminidase 3, membrane sialidase

NCBI Official Synonym Symbols

SIAL3

NCBI Protein Information

sialidase-3

UniProt Gene Name

NEU3

UniProt Protein Name

Sialidase-3

UniProt Synonym Protein Names

Ganglioside sialidase; Membrane sialidase; N-acetyl-alpha-neuraminidase 3

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Sialidase-3 (NEU3), ELISA Kit

Cat PEK05299

UniProt Entry Name

NEUR3_HUMAN

UniProt Secondary Accession

Q9NQE1; A8K327

UniProt Related Accession

Q9UQ49

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY