

Hyaluronan Synthase 3 (HAS3), ELISA Kit

Cat PEK05477

Species

Plant

Full Product Name

Plant Hyaluronan Synthase 3 (HAS3) ELISA Kit

Product Gene Name

HAS3 elisa kit

Molecular Weight

31,200 Da

Storage

Store all reagents at 2-8°C

Protein Family

Hyaluronan synthase

NCBI Accession

AAK73797.1

NCBI GI

14915764

NCBI GeneID

3038

NCBI Official Full Name

hyaluronan synthase 3

NCBI Official Symbol

HAS3

NCBI Official Synonym Full Names

hyaluronan synthase 3

NCBI Protein Information

hyaluronan synthase 3

UniProt Gene Name

HAS3

UniProt Synonym Gene Names

HA synthase 3

UniProt Protein Name

Hyaluronan synthase 3

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Hyaluronan Synthase 3 (HAS3), ELISA Kit

Cat PEK05477

Hyaluronate synthase 3; Hyaluronic acid synthase 3; HA synthase 3

UniProt Entry Name

HYAS3_HUMAN

UniProt Secondary Accession

Q8WTZ0; Q9NYP0; A8K5T5

UniProt Related Accession

O00219

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY