

Late Cornified Envelope Protein 3E (LCE3E), ELISA Kit

Cat PEK05582

Species

Plant

Full Product Name

Plant Late Cornified Envelope Protein 3E (LCE3E) ELISA Kit

Product Gene Name

LCE3E elisa kit

Molecular Weight

9,507 Da

Storage

Store all reagents at 2-8°C

Protein Family

Late cornified envelope protein

NCBI Accession

NP_848522.1

NCBI GI

30410033

NCBI GeneID

353145

NCBI Official Full Name

late cornified envelope protein 3E

NCBI Official Symbol

LCE3E

NCBI Official Synonym Full Names

late cornified envelope 3E

NCBI Official Synonym Symbols

LEP17

NCBI Protein Information

late cornified envelope protein 3E

UniProt Gene Name

LCE3E

UniProt Synonym Gene Names

LEP17

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Late Cornified Envelope Protein 3E (LCE3E), ELISA Kit

Cat PEK05582



Late cornified envelope protein 3E

UniProt Synonym Protein Names

Late envelope protein 17

UniProt Entry Name

LCE3E_HUMAN

UniProt Secondary Accession

A2RRM6

UniProt Related Accession

Q5T5B0

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
