

Plasminogen Activator (PLGA), ELISA Kit

Cat PEK05587

Species

Plant

Full Product Name

Plant Plasminogen Activator (PLGA) ELISA Kit

Product Gene Name

PLGA elisa kit

Molecular Weight

44,373 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

AAA98809.1

NCBI GI

339818

NCBI GenelD

5327

NCBI Official Full Name

plasminogen activator

NCBI Official Symbol

PLAT

NCBI Official Synonym Full Names

plasminogen activator, tissue type

NCBI Official Synonym Symbols

TPA; T-PA

NCBI Protein Information

tissue-type plasminogen activator

UniProt Gene Name

PLAT

UniProt Synonym Gene Names

t-PA; t-plasminogen activator; tPA

UniProt Protein Name

Tissue-type plasminogen activator

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Plasminogen Activator (PLGA), ELISA Kit

Cat PEK05587

INN: Alteplase; INN: Reteplase

UniProt Entry Name

TPA_HUMAN

UniProt Secondary Accession

Q15103; Q503B0; Q6PJA5; Q7Z7N2; Q86YK8; Q9BU99; Q9BZW1; A8K022; B2R8E8

UniProt Related Accession

P00750

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY